

The Hall® UltraPower® Surgical Drill System



endocon[®]

The UltraPower® Surgical Drill System

your first choice for lasting power

The UltraPower Surgical Drill System is the ultimate solution for addressing a wide variety of demanding procedures and the need for a more economical and convenient OR strategy.

Optimum site management, with built-in economy

- High Speed Drill has a max. speed of 70,000 rpm and 7.5 in. oz. of stall torque
- Low Speed Drill has a max. speed of 20,000 rpm and 23 in. oz. of stall torque
- Provides precision cuts quickly and cleanly without bogging down
- Perforator Drive handpiece accepts any cranial perforator with a Hudson shank
- Irrigation System available
- Convenient variable-speed foot control ensures control at surgical site by allowing you to choose the speed that makes the most sense for the situation
- Operating at an economical 100 psi, the lubeless UltraPower System is dependable during even the longest of procedures

The UltraPower Perforator Drive



The UltraPower High Speed Drill

The UltraPower Low Speed Drill

Connect with superior cutter technology

- Wide array of burs
- Diversely sized and shaped heads with assorted shaft lengths ensure an ideal choice for each procedure
- Positive lock collet mechanically secures the bur for extra convenience and safety
- Wide variety of angled attachments and bur guards accommodate various surgical requirements
- With over 30 years of experience in pneumatics, Hall Surgical is proud to show you the essence of choice, putting the power of versatility at your fingertips

The UltraPower High Speed Drill with Universal Hose and 20° Angle Attachment



The UltraPower Foot Control

Complete selection of Bur Guards and Angled Attachments





endocon[©]

endocon GmbH
Carl-Bosch-Straße 4
69115 Heidelberg

Simplify Mobility

T +49 6221 7575 55
F +49 6221 7575 45
E info@endocon.de
I www.endocon.de